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

Figure 1

The structures of levofloxacin (A) and amlodipine besylate (B).


Figure 2


Typical calibration plot in the range 300–10000 ng mL<sup>-1</sup>  
( $Y = 6.132 + 0.015X$ ;  $R = 0.998$ ,  $sdv = 5.15$ ).

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

▲ Example of a HPTLC plate with the violet amino-propanol zone; respective sample and standard tracks were developed from both plate sides in the HDC and just heated for derivatization

▲ 3D-graphic of the close-up range of the amino-propanol peak on different tracks measured at 486 nm (five-point calibration is marked)

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
Validation parameter	Results
LOD	4.5 µg/mL
LOQ	15 µg/mL
<i>(both related to the applied volume 2 µL)</i>	
Linearity (coefficient of correlation)	0.9979
Mean recovery	102 %
Mean repeatability (RSD)	±4.9 %
<i>(performed at 5 concentration levels with n= 3 for each calibration level)</i>	
Intermediate precision (RSD, n= 9)	±5.7 %

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
Calibration for acetylsalicylic acid, 150-350 ng (peak height)

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

Track 2 , ID: S7897-01 no 1 capsule in 50 ml Tol diluted 1:40

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
X(Average)			Target content T (label claim)	30 ng
29.53 ng			Average (calculated)	30.08 ng
			Mean of individual contents expressed as a percentage of the label claim $\bar{X}$ (calcul.)	$\bar{X} = 100.25$
X(Average)	X(Average)	X(Average)	Sample standard deviation s (calcul.)	s = 2.76
31.28 ng	30.43 ng	28.64 ng	Acceptability constant k (n=10)	k = 2.4
			M (case 1), if $98.5\% \leq \bar{X} \leq 101.5\%$	M = $\bar{X}$
X(Average)	X(Average)	X(Average)	Acceptance value AV	AV = $ M - \bar{X}  + ks$
30.17 ng	30.99 ng	30.35 ng	AV (calcul.)	<b>AV = 6.6</b>
			Max. allowed AV, L1 = 15.0	15.0
X(Average)	X(Average)	X(Average)		
30.05 ng	30.35 ng	28.97 ng		

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


Image white light


Image white light (enhanced)

- 1: Bilberry dry extract
- 2: Bilberry dry extract spiked with amaranth (spiking level 0.25 %)
- 3: Amaranth

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1: Arachis oil  
2: Sesame oil  
3: Sesame oil roasted  
4: Maize oil  
5: Linseed oil

6: Sunflower oil  
7: Almond oil  
8: Wheat germ oil  
9: Evening primrose oil  
10: Safflower oil type 1

11: Soy bean oil  
12: Borage oil  
13: Argan oil  
14: Grape seed oil

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
Height calibration for substance: Hydrastine

POLYNOMIAL REGRESSION Calibration method

This calibration is based on 7 standards.

Calibration equation (Y = height X = amount detected (ng))

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


Standards


Cholesteryl oleate  
Palmitic acid /  
Tetracosanoic acid  
Cholesterol  
Ceramide NS  
Ceramide NP  
Ceramide AP  
Cholesteryl-3-sulfate

Detection:  
CuSO<sub>4</sub> in 8% phosphoric  
acid, 150°, 20min

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
Densitogram of eight phospholipids

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

Calibration function for phosphatidyl-choline (PC)

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
Testdye III, application volume 1, 2 and 3 uL

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

- High throughput (46 runs in 15 min by (anti-)parallel development, 15 min-staggered offline system) → 1000 samples/8h-day
- Resulting in 20 s per sample with 330  $\mu$ L solvent consumption
- Almost no disposal costs < 0.01 Cent/sample
- Selective derivatization → compensates low separation power
- Reduced sample preparation: no SPE
- Analysis without acetonitrile!

- Ultra-rapid HPLC with 2 min gradient: 720 runs/24-h day
- Sample preparation: Need of SPE for MS or ELSD as detector

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
The screenshot shows the 'winCATS Image Comparison' software window. At the top, there is a menu bar with 'File' and 'Help'. Below it is a toolbar with icons for back, forward, and search. The main area contains a table with the following columns: 'Show', 'Image', 'Refer...', 'Sample ID', 'Track', and 'Source ID'. The table lists several samples, including 'Rheum officinale', 'Rhei sinensis', 'Rheum rhaponticum', and 'Rheum palmatum'. Below the table, there is a grid of gel images. The first three images (labeled 19, 15, and 3) are highlighted in blue, indicating they are used for comparison. The rest of the images are labeled with their respective sample IDs. At the bottom of the window, there is a status bar with fields for 'C:\CAMAG\winCATS\Data\ Modified', 'Zoom: 1.00', '0713W006', and 'Jan Masthoff'.

Show	Image	Refer...	Sample ID	Track	Source ID
<input checked="" type="checkbox"/>	0	<input type="checkbox"/>	Rheum officinale	1	C:\CAMAG\winC...
<input checked="" type="checkbox"/>	9	<input type="checkbox"/>	Rhei sinensis	2	C:\CAMAG\winC...
<input checked="" type="checkbox"/>	17	<input type="checkbox"/>	Rhei sinensis	2	C:\CAMAG\winC...
<input checked="" type="checkbox"/>	1	<input type="checkbox"/>	Rhei sinensis	2	C:\CAMAG\winC...
<input checked="" type="checkbox"/>	11	<input type="checkbox"/>	Rheum rhaponticum	4	C:\CAMAG\winC...
<input checked="" type="checkbox"/>	2	<input type="checkbox"/>	Rheum palmatum	3	C:\CAMAG\winC...
<input checked="" type="checkbox"/>	3	<input checked="" type="checkbox"/> isRef	Rheum rhaponticum	4	C:\CAMAG\winC...
<input checked="" type="checkbox"/>	4	<input type="checkbox"/>	Rheum officinale	5	C:\CAMAG\winC...

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CAMAG BIBLIOGRAPHY SERVICE PLAVAR CHROMATOGRAPHY CBS 101 • September 2008

Automatic Sample Application  
1974 and today

**CBS**  
CAMAG Bibliography Service

Excerpts from CBS 51 - 99 • Keyword:

**Pharmaceutical and biomedical applications**

J. KRZEK\*, U. HUBICKA, J. SZCZEPANCZYK (\*Jagiellonian University, Collegium Medicum, Department of Inorganic and Analytical Chemistry, Medyczna 9, 30-688 Krakow, Poland): **High-performance thin-layer chromatography with densitometry for the determination of ciprofloxacin and impurities in drugs.** J. Assoc. Off. Anal. Chem. 88, 1530-1535 (2005). HPTLC of ciprofloxacin and degradation products (an ethylenediamine compound, a desfluoro compound, a by-compound, and fluoroquinolonic acid) on silica gel with chloroform - methanol - 25 % ammonia 43:43:14. Quantitative determination by densitometric analysis at 330 nm for fluoroquinolonic acid and at 277 nm for the other compounds. The method showed high sensitivity (limit of detection 10 to 44 ng), a wide linearity range (3 to 20 µg/mL), and good precision (2.32 to 6.46 % relative standard deviation) and accuracy (recovery rates 98.6 to 101.5 %) for individual constituents.

quality control, densitometry, quantitative analysis, HPTLC

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INSTRUMENTS AND TOOLS FOR ALL STEPS  
OF MODERN THIN-LAYER CHROMATOGRAPHY

[My Account](#) | [Search](#)

[Products](#)

[Support](#)

[News & Events](#)

[Literature](#)

[About CAMAG](#)

[Laboratory](#)

**What's new?**

**New CAMAG TLC-MS Interface:**  
For direct extraction of compounds from TLC/HPTLC layers into MS.

[More information >>](#)

[Press release >>](#)

**From the Lab**

**New publications**  
Validated HPTLC method for the identification of Hoodia gordonii >>

Validation of HPTLC methods for the identification of botanicals in a cGMP environment >>

**Did you know?**

CAMAG manufactures **Aluminium oxide** which is one of the most versatile sorbents for preparative chromatography. Various types of Alox are used in

- Preparative Column Chromatography
- Flash Chromatography
- Selective Filtration

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my Account | Search

Services
Methods
Publications
Virtual Classroom
About us
CAMAG home

**Methods**

A broad range of

List of all applicati

List of all applicati

New

Application Notes

Validated Methods

Medicinal Plants

Food

Pharma / Biotechnology / Environmental / Forensic

Registered users can download the Application Notes in pdf-format. Choose by subject or use search.

- Medicinal Plants
- Food
- Pharma / Biotechnology / Environmental / Forensic / Cosmetics

**New**

ID	Application Note
A-95.1	<p><b>Determination of metronidazole in pharmaceutical dosage form by HPTLC</b></p> <p>This method is suitable for the quantification of metronidazole in pharmaceutical dosage form. Metronidazole is a nitroimidazole antibiotic marketed as a prescription drug in various countries.</p>
A-91.1	<p><b>Determination of phospholipids by HPTLC</b></p> <p>Phospholipids are important compounds in biological systems and of high interest to the pharmaceutical and cosmetic industry. Quantitative determination of phospholipids is also important for the food industry with respect to declaration of prophylactic and therapeutic effects of dietary phospholipids. This method is suitable for screening and quantification of lysophosphatidic acid, lysophosphatidylcholine, lysophosphatidylethanolamine, lysophosphatidylinositol, phosphatidic acid, phosphatidyl choline, phosphatidyl ethanolamine and phosphatidylinositol</p>

APPLICATION NOTES

**HPTLC Identification of Common Horse Tail**  
*(Equisetum arvense)* F-12.1

**Scope:** This HPTLC method identifies *Equisetum arvense* (Common Horse Tail) and allows discrimination from the adjacent *Equisetum sylvaticum* at the 25% level. The method is based on a collaborative study of Official Medicines Control Laboratories (OMCL) and features fast and simplified sample preparation.

**Sample:** 1.0 g of the dried drug is mixed with 10 ml of methanol, sonicated for 10 min and centrifuged. The supernatant is used as test solution.

**Standards (optional):** 5 mg of each, rutin, hyperoside, and quercetin, and are individually dissolved in 10 mL of methanol.

**Default reagents:** Natural Products reagent (NP reagent) 1 g of dibutyltinous acid and acetylacetone dissolved in 200 mL of ethyl acetate.  
Microgel reagent: 10 g of polyethylene glycol 400 (Carbowax) are dissolved in 200 mL of dichloromethane.

**Chromatographic conditions:**

**Stationary phase:** HPTLC plates silica gel 60 F254 (Merck), 10x10 cm or 20x10 cm.

**Mobility Phase:** Formic acid, acetic acid, water, ethyl acetate (15:15:20:134).

**Sample application:** 10 µL of test solution and 2 µL of standards are applied as 8 mm bands, min. 2 mm apart, 8 mm from lower edge of plate.

**Development:** 10x10 cm or 20x10 cm Thin Layer Chamber, saturated for 20 min. Other papers, 6 mL (depending 10 mL) developing solvent per trough, developing distance 70 mm from lower

Please contact CAMAG for more application notes and products for analysis of natural products  
www.camag-laboratory.com 1 of 2 11/07

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### **Validation of High-Performance Thin-Layer Chromatographic Methods for the Identification of Botanicals in a cGMP Environment**

EIKE REICH and ANNE SCHIBLI

CAMAG Laboratory, Sonnenmattstrasse 11, CH-4132 Muttenz, Switzerland

ALISON DeBATT

CAMAG Scientific Inc., 515 Cornelius Hamett Dr, Wilmington, NC 28401

---

**Current Good Manufacturing Practices (cGMP) for botanicals stipulates the use of appropriate methods for identification of raw materials. Due to natural variability, chemical analysis of plant material is a great challenge and requires special approaches. This paper presents a comprehensive proposal to the process of validating qualitative**

system to ensure safety, quality, and consistency of products. Because quality can not be tested into a product, it is important that quality is already an attribute of the raw material of any given production process. For “botanicals,” which can be fresh, or processed plants, extracts, or other preparations, the identity is one of the central elements of quality. cGMP

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Sonnenmattstrasse 11  
CH-4132 Muttenz  
Lab@camag.com

### Request for Service

<b>Information about customer</b>	
Name of company:	
Person to contact:	
Department:	
Address:	
Phone:	
Fax:	
E-mail:	
Field of activity/business:	
<b>Information about TLC problem</b>	
Evaluation:	<input type="checkbox"/> qualitative <input type="checkbox"/> quantitative
Type of service requested:	<input type="checkbox"/> literature research only <input type="checkbox"/> feasibility study <input type="checkbox"/> method development <input type="checkbox"/> validation
Current method of analysis (please attach method):	
<input type="checkbox"/> TLC/HPTLC <input type="checkbox"/> HPLC <input type="checkbox"/> GC <input type="checkbox"/> other <input type="checkbox"/> none	
Type of sample:	
Matrix:	<input type="checkbox"/> pharmaceutical formulation (tablet, syrup, ) <input type="checkbox"/> cosmetic <input type="checkbox"/> plant extract <input type="checkbox"/> oil / glycerin

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